


FlexTip Bougie

Traditional bougies are valuable instruments, especially for difficult cases; but sometimes they're just not enough. Bougie tip control makes navigating the airway much easier.

Whether you use DL or VL, the FlexTip Bougie will help find the way. Simply use your thumb to move the tabs of the integrated slider up and down to flex or retro-flex the tip of this innovative new bougie from P3 Medical.

- Flexible, steerable tip makes it easier to navigate difficult airways
- Tabs are easy to grip and manoeuvre even with gloves
- Slider tabs situated so they will not need to pass through vocal cords
- Blunt silicone tip reduces the risk of tracheal trauma
- Bright phosphorous coating on tip makes it highly visible even with the ultra-violet used by some laryngoscopes
- Introduce with anterior flexion, straighten to advance the tube or retro-flex if needed
- Graduation marks from 10 to 50cm, 65cm overall
- Improves intubation, saving time, leading to fewer complications
- Latex free technology



Product Code: **CM001-EU**

Individually packed sterile, 10 units per box

Patents pending

Svensk distributör:

 **Medu**

E-post: info@medu.se
Tel: 08-31 00 46
www.medu.se

P3 Medical Limited

October 2017